


Minimize your balloon profile to maximize your treatment possibilities.



Pliaform™ balloon-surface texture

Advance Micro™ 14

ULTRA LOW-PROFILE PTA BALLOON CATHETER

Microtherapy

The ultra low profile of these balloons allows you to attempt interventions that you never thought were possible.

Ideal for Below the Knee

These balloons were designed specifically to treat vessels from the knee all the way to the digital arteries.


Reduced Friction

The Pliaform balloon-surface texture* and hydrophilic coating reduce friction to help you cross lesions.

*on all sizes except 1.5 mm diameters



www.cookmedical.com



Actual size

Advance Micro™ 14

ULTRA LOW-PROFILE PTA BALLOON CATHETER

The Advance Micro 14 is used for percutaneous transluminal angioplasty (PTA) of lesions in peripheral arteries, including iliac, renal, popliteal, posterior tibial, anterior tibial, peroneal, pedal, femoral, iliofemoral, radial, brachial, ulnar and internal pudendal arteries, as well as of obstructive lesions of native or synthetic arteriovenous dialysis fistulas.

- Ultra-low profiles allow these balloons to pass through a 2.9 Fr ID Micropuncture® Pedal Access Set introducer.
- All of the balloon sizes are compatible with 3 Fr sheaths.
- The hydrophilic coating enhances pushability and trackability.
- Balloon lengths of up to 12 cm allow the treatment of longer lesions with fewer inflations and deflations.
- The nominal pressure for all sizes is 8 atm. The rated burst pressure is 16 atm.

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
50 cm Shaft – .014 inch Wire Guide Compatible			
G26090	PTA3-14-50-1.5-2	1.5	2
G26091	PTA3-14-50-1.5-4	1.5	4
G26092	PTA3-14-50-2-2	2	2
G26093	PTA3-14-50-2-3	2	3
G26094	PTA3-14-50-2-4	2	4
G26095	PTA3-14-50-2-6	2	6
G26096	PTA3-14-50-2-8	2	8
G26098	PTA3-14-50-2-10	2	10
G26099	PTA3-14-50-2-12	2	12
G26100	PTA3-14-50-2.5-2	2.5	2
G26101	PTA3-14-50-2.5-3	2.5	3

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
G26102	PTA3-14-50-2.5-4	2.5	4
G26103	PTA3-14-50-2.5-6	2.5	6
G26104	PTA3-14-50-2.5-8	2.5	8
G26105	PTA3-14-50-2.5-10	2.5	10
G26106	PTA3-14-50-2.5-12	2.5	12
G26107	PTA3-14-50-3-2	3	2
G26109	PTA3-14-50-3-3	3	3
G26110	PTA3-14-50-3-4	3	4
G26111	PTA3-14-50-3-6	3	6
G26112	PTA3-14-50-3-8	3	8
G26113	PTA3-14-50-3-10	3	10
G26115	PTA3-14-50-3-12	3	12

Reference Part Number Key

PTA3-14-50-1.5-2



Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
90 cm Shaft – .014 inch Wire Guide Compatible			
G26117	PTA3-14-90-1.5-2	1.5	2
G26118	PTA3-14-90-1.5-4	1.5	4
G26119	PTA3-14-90-2-2	2	2
G26121	PTA3-14-90-2-3	2	3
G26122	PTA3-14-90-2-4	2	4
G26123	PTA3-14-90-2-6	2	6
G26195	PTA3-14-90-2-8	2	8
G26400	PTA3-14-90-2-10	2	10
G26401	PTA3-14-90-2-12	2	12
G26402	PTA3-14-90-2.5-2	2.5	2
G26403	PTA3-14-90-2.5-3	2.5	3
G26404	PTA3-14-90-2.5-4	2.5	4
G26405	PTA3-14-90-2.5-6	2.5	6
G26406	PTA3-14-90-2.5-8	2.5	8
G26407	PTA3-14-90-2.5-10	2.5	10
G26601	PTA3-14-90-2.5-12	2.5	12
G26602	PTA3-14-90-3-2	3	2
G26604	PTA3-14-90-3-3	3	3
G26605	PTA3-14-90-3-4	3	4
G26606	PTA3-14-90-3-6	3	6
G26607	PTA3-14-90-3-8	3	8
G26608	PTA3-14-90-3-10	3	10
G26610	PTA3-14-90-3-12	3	12

Order Number	Reference Part Number	Inflated Balloon Diameter mm	Balloon Length cm
150 cm Shaft – .014 inch Wire Guide Compatible			
G26611	PTA3-14-150-1.5-2	1.5	2
G26614	PTA3-14-150-1.5-4	1.5	4
G26620	PTA3-14-150-2-2	2	2
G26621	PTA3-14-150-2-3	2	3
G26622	PTA3-14-150-2-4	2	4
G26625	PTA3-14-150-2-6	2	6
G26626	PTA3-14-150-2-8	2	8
G26636	PTA3-14-150-2-10	2	10
G26641	PTA3-14-150-2-12	2	12
G26643	PTA3-14-150-2.5-2	2.5	2
G26646	PTA3-14-150-2.5-3	2.5	3
G26653	PTA3-14-150-2.5-4	2.5	4
G26665	PTA3-14-150-2.5-6	2.5	6
G26677	PTA3-14-150-2.5-8	2.5	8
G26678	PTA3-14-150-2.5-10	2.5	10
G26679	PTA3-14-150-2.5-12	2.5	12
G26680	PTA3-14-150-3-2	3	2
G26681	PTA3-14-150-3-3	3	3
G26706	PTA3-14-150-3-4	3	4
G26707	PTA3-14-150-3-6	3	6
G26710	PTA3-14-150-3-8	3	8
G26719	PTA3-14-150-3-10	3	10
G26724	PTA3-14-150-3-12	3	12

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Compatible products


ACCESS

The **Micropuncture® Pedal Access Set** can facilitate a retrograde tibiopedal approach, which is a new revascularization technique that aims to give you another chance at limb salvage.*

TARGET

When confronted with a chronic total occlusion in the legs, you can harness controlled power by choosing from four tip-load options of the **Approach® CTO Microwire Guide**.

The hydrophilic tip and supportive shaft of the **Approach® Hydro ST Microwire Guide** can help you target anatomy that is tortuous or below the knee.



Approach® CTO
Microwire Guide



Approach® Hydro ST
Microwire Guide




Micropuncture®
Pedal Access Set

The **Shuttle® Tibial Guiding Sheath** pairs with specially designed **Slip-Cath® Angiographic Catheters** to simplify access while minimizing device exchanges.

The **CXI® Support Catheter** provides the strength, visibility, lubricity and low profile that you need to help you reach and cross lesions.



Shuttle® Tibial
Infrapopliteal
Access System



CXI®
Support Catheter

* Montero-Baker M, Schmidt A, Bräunlich S, et al. Retrograde approach for complex popliteal and tibioperoneal occlusions. *J Endovasc Ther.* 2008;15(5):594-604.

Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
E-mail: cau.custserv@cookmedical.com

