



Cross lesions with just one set



Four configurations and
two lengths available



TriForce[®] PERIPHERAL CROSSING SET

Optimized pushability


Designed coaxial pairing of CXI[®] catheter with Flexor[®] sheath facilitates advancement and maneuverability.

Enhanced visibility

Distal ends of catheter and sheath have increased radiopacity for visualization.

Trackability

Smooth transition of coaxial system combined with hydrophilic coating supports crossing the lesion.



TriForce®

PERIPHERAL CROSSING SET

- TriForce Peripheral Crossing Set includes an outer 5.0 Fr Flexor® sheath, a 4.0 Fr CXI® Support Catheter, and a Peel-Away® sheath.
- Flexor outer sheath tapers to a 4.0 Fr inner diameter at the tip to ensure a smooth transition to the CXI Support Catheter.
- Flexor outer sheath features radiopaque tip and hydrophilic coating.
- Accepts a .035 inch diameter wire guide.

Order Number	Reference Part Number	Catheter Fr/Length cm	Sheath ID/Length Fr/cm	Sheath Tip Fr	Accepts Wire Guide Diameter Inch
Set with Straight Tip Sheath and Straight Tip Catheter					
G56412	KCXS-5.0-35-65-RB-0/0-HC	4.0/65	5.0/55	4.0	.035
G56416	KCXS-5.0-35-100-RB-0/0-HC	4.0/100	5.0/90	4.0	.035
Set with Straight Tip Sheath and Curved Tip Catheter					
G56413	KCXS-5.0-35-65-RB-0/DAV-HC	4.0/65	5.0/55	4.0	.035
G56417	KCXS-5.0-35-100-RB-0/DAV-HC	4.0/100	5.0/90	4.0	.035
Set with Curved Tip Sheath and Straight Tip Catheter					
G56414	KCXS-5.0-35-65-RB-MPB/0-HC	4.0/65	5.0/55	4.0	.035
G56418	KCXS-5.0-35-100-RB-MPB/0-HC	4.0/100	5.0/90	4.0	.035
Set with Curved Tip Sheath and Curved Tip Catheter					
G56415	KCXS-5.0-35-65-RB-MPB/DAV-HC	4.0/65	5.0/55	4.0	.035
G56419	KCXS-5.0-35-100-RB-MPB/DAV-HC	4.0/100	5.0/90	4.0	.035

Set includes 5.0 Fr Flexor outer sheath, 4.0 Fr CXI Support Catheter and Peel-Away sheath.

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 8007615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: customersupport@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-8.5X11