











Catheters | Features

		PTFE/Polyethylene guide catheter	Polyethylene guide catheter	PTFE guide catheter	MicroJol™ technology	Support cannula	Distal ink marks	Proximal ink marks	EchoTip®	Stiff obturator available	Malleable obturator available	Bulb tip	Precurved
Guardia™ Access Embryo Transfer Catheter G34782 K-JETS-6019 G34783 K-JETS-7019			•		•	•		•		•	•	•	•
Guardia™ AccessET Curved Embryo Transfer Catheter G53247 K-JETS-6019-ET G53007 K-JETS-7019-ET			•		•	•		•	•	•	•	•	•
Guardia™ Access Embryo Transfer Catheter with Internal Support Cannula G02346 K-JETS-7019-INT			•		•	•		•		•	•	•	•
Guardia™ Access Nano Embryo Transfer Catheter G24216 K-JETS-551910-S		•		•	•	•	•	•			•	•	•
Pivot Embryo Transfer Catheter Set G30457 K-PETS-5000					•		•	•		•	•		














Other embryo transfer catheter products

Echosight® Jansen-Anderson Intratubal Transfer Set G17946 K-J-JITS-572900	
Towako Transmyometrial Embryo Transfer Set G18742 K-TTET-19-32.5-TOWAKO	
Guardia™ Access Malleable Obturator G34784 K-JET-7002	
Guardia™ Obturator G54189 J-UOB-2828 Malleable G54190 J-UOB-4028 Malleable G54191 J-UOB-4028-ST Stiff	
Soft-Pass™ Malleable Obturator G53744 J-SP-4420	

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Choosing embryo transfer catheters

Catheters | Features

		PTFE/Polyethylene guide catheter	Polyethylene guide catheter	PTFE guide catheter	Microvol [™] technology	Support cannula	Distal ink marks	Proximal ink marks	EchoTip [®]	Stiff obturator available	Malleable obturator available	Bulb tip	Precurved
Soft-Pass™ Embryo Transfer Catheter Set G17934 K-J-SPPE-681710 with support G17932 K-J-SPPE-681200 without support G17933 K-J-SPPE-681700 without support			•			•	•	•		•	•		
			•				•	•		•	•		
			•				•	•		•	•		
Soft-Pass™ Embryo Transfer Catheter Set with EchoTip[®] G18954 K-J-SPPE-681210-ET with support G18955 K-J-SPPE-681710-ET with support G18657 K-J-SPPE-681200-ET without support G18658 K-J-SPPE-681700-ET without support			•			•	•	•	•	•	•		
			•			•	•	•	•	•	•		
			•				•	•	•	•	•		
			•				•	•	•	•	•		
Soft-Pass™ Mock Transfer Catheter Set G30399 J-SPPE-681710-ET-MC G17678 J-SPPE-681200-MC G17679 J-SPPE-681700-MC			•			•	•	•	•	•	•		
			•				•	•		•	•		
			•				•	•		•	•		
Soft-Pass™ Embryo Transfer Catheter G17931 K-J-SP-681710 with support G17902 K-J-SP-681700 without support				•		•	•	•		•	•		
				•			•	•		•	•		
Soft-Pass™ Embryo Transfer Catheter with EchoTip[®] G27206 K-J-SP-681700-ET				•			•	•	•	•	•		

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4