

Place embryos
with confidence.



Guardia™
CURVED EMBRYO TRANSFER CATHETERS



Guardia™

CURVED EMBRYO TRANSFER CATHETERS

Guardia embryo transfer catheters are designed for ease of access:

- The bulb-shaped tip promotes passage through the cervix and into the uterine cavity.
- The curved guide catheter and soft, flexible transfer catheter may aid in navigation.
- The Microvol™ technology decreases the volume of media that is required for transfer.





Guardia™ Access

CURVED EMBRYO TRANSFER CATHETER

Cook Medical's original curved embryo transfer catheter.

The original



Guardia™ AccessET

CURVED EMBRYO TRANSFER CATHETER

The EchoTip echogenic tip enhances visibility under ultrasound.

With EchoTip™


Smaller diameter



Guardia™ Access Nano

CURVED EMBRYO TRANSFER CATHETER


The smaller diameter provides navigability and flexibility.



Guardia™

CURVED EMBRYO TRANSFER CATHETERS

Guardia curved embryo transfer catheters are used to place in vitro fertilized (IVF) embryos into the uterine cavity.



Order Number	Reference Part Number	Guide Catheter Fr	Guide Catheter Length cm	Transfer Catheter Fr	Transfer Catheter Length cm	Guide Catheter Hub Type
Guardia Access						
Catheters Preassembled						
G34782	K-JETS-6019	6.6	16.7	2.8	23	small
Catheters Packaged Separately						
G34783	K-JETS-7019	6.6	17.3	2.8	24	large
Guardia AccessET						
G53247	K-JETS-6019-ET	6.6	16.7	2.8	23	small
G53007	K-JETS-7019-ET	6.6	17.3	2.8	24	large
Guardia Access Nano						
G24216	K-JETS-551910-S	5.5	17.3	2.8	24	large
Optional Malleable Obturator						
G34784	K-JET-7002					
Optional Malleable Adjustable Obturator						
G54189	J-UOB-2828					

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4