



Selectively catheterize and cannulate the proximal fallopian tube.



The transparent catheter tip allows visualization of the inner catheter and wire guide.

Modified Novy

CORNUAL CANNULATION SET




Modified Novy

CORNUAL CANNULATION SET

Intended for use through the operating channel of a hysteroscope or other uterine access device, for hysteroscopic or fluoroscopic selective catheterization and cannulation of the proximal fallopian tube(s), followed by the introduction of chromotubation solution or contrast medium, in the evaluation of tubal patency.

The wire guide in the Modified Novy Cornual Cannulation Set is intended only to facilitate placement of the inner catheter. It is not intended for tubal recanalization and should not be advanced beyond the tubal isthmus.

- The Torcon® polyethylene introducing catheter facilitates access to the cornual os.
- The transparent catheter tip allows visualization of the inner catheter and wire guide.



| Order Number | Reference Part Number | Introducing Catheter Fr | Introducing Catheter Length cm | Inner Catheter Fr | Inner Catheter Length cm | Wire Guide Diameter inch/mm |
|--------------|-----------------------|-------------------------|--------------------------------|-------------------|--------------------------|-----------------------------|
| G17478 | J-NCS-503570 | 5.0 | 35 | 3.0 | 50 | .018/.46 |
| G17558 | J-NCS-504070 | 5.0 | 40 | 3.0 | 55 | .018/.46 |

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4