


For urinary tract decompression in fetuses 18-32 weeks.*



*For the treatment of fetal postvesicular obstructive uropathy.

Illustration by Lisa Clark


HARRISON FETAL BLADDER STENT SET



HARRISON FETAL BLADDER STENT SET

Used in fetal urinary tract decompression following the diagnosis of fetal postvesicular obstructive uropathy in fetuses of 18 to 32 weeks gestational age.

- The EchoTip® technology enhances the visualization of the needle tip under ultrasound.
- The usable length between pigtails is 1.5-3.5 cm.



13 gage, 18 cm long EchoTip trocar needle



.038 inch (.97 mm) diameter, 40 cm long PTFE-coated stainless steel wire guide



5.0 Fr, 24 cm long positioner



Multilength double-pigtail stent

Order Number	Reference Part Number	Guide Needle gage	Guide Needle Length cm	Wire Guide Diameter inch/mm	Wire Guide Length cm	Multilength Stent Fr
G16766	J-HFBS-503540	13	18	.038/.97	40	5.0

Humanitarian device: Authorized by federal (USA) law for use in the treatment of fetal postvesicular obstructive uropathy. The effectiveness of this device for this use has not been demonstrated.

CAUTION: Federal (USA) law restricts this device to sale, distribution, and use by or on the order of a physician with appropriate training and experience.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
 EDI: cookmedical.eu/edi
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001702, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Iceland: +354 800 7615, IS.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com
 EDI: cookmedical.com/edi.do
Americas:
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com
Australia:
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
 E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4