







Cook Medical offers a wide range of catheters with placement and access features that support your practice and preference during procedures.

Catheters | Features



		Balloon	Acorn-shaped positioner	Working channel	Removable stylet	Support sheath	Formable catheter	Sonohysterography	Hysterosalpingography	Selective salpingography	Endometrial biopsy	Intrauterine insemination
Soules Intrauterine Catheter G16473 J-IUIO-532500 5.3 Fr, 25 cm long								•				•
Goldstein Sonohysterography Catheter G17149 J-GSHC-532600 5.3 Fr, 26 cm long G19012 J-GSHC-542600-SV 5.4 Fr, 26 cm long			•				•	•				
Goldstein SonoBiopsy™ Catheter G46190 J-GBSX-072026 7.2 Fr, 26 cm long			•				•	•			•	
Margolin HSG Cannula G16504 J-MHSG-442500 16.0 Fr, 25 cm long			•		•				•			
Cook® Silicone Balloon HSG Catheter G17147 J-CHSG-503000 5.0 Fr, 30 cm long G17870 J-CHSG-703000 7.0 Fr, 30 cm long		•					•	•	•			
Intrauterine Access Balloon Catheter G30397 J-IAB-090023 9.0 Fr, 23 cm long		•		•	•				•	•		

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Choosing diagnostic catheters

Cook Medical offers a wide range of catheters with placement and access features that support your practice and preference during procedures.

Catheters | Features

		Balloon	Acorn-shaped positioner	Working channel	Removable stylet	Support sheath	Formable catheter	Sonohysterography	Hysterosalpingography	Selective salpingography	Endometrial biopsy	Intrauterine insemination
Selective Salpingography Catheter G17359 J-SSG-554086 5.5 Fr, 40 cm long Tip angle: 65 degrees				•						•		
Rösch-Thurmond Fallopian Tube Catheterization Set G04324 FTC-900 5.5 Fr torque control catheter, 3.0 Fr injection catheter, .015 inch diameter Cope mandril wire guide				•					•	•		

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

Customer Service

EU Website: cookmedical.eu
EDI: cookmedical.eu/edi
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001702, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Iceland: +354 800 7615, IS.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Norway: +47 23162968, no.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

USA Website: cookmedical.com

EDI: cookmedical.com/edi.do

Americas:

Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335

E-mail: orders@cookmedical.com

Australia:

Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283

E-mail: cau.custserv@cookmedical.com



AI, ESC, IR, OHNS, PI, RH, SUR-A4