



Amplatz

FIXED CORE WIRE GUIDE




Access and navigate the ureter.



Amplatz

FIXED CORE WIRE GUIDE



Amplatz wire guides are used to gain ureteral access, to establish a tract, and to assist in the placement, replacement and exchange of medical devices during urological procedures.

Intended for one-time use. These wire guides are not intended for PTCA use.

Global Product Number	Order Number	Diameter inch	Length cm	Tip Configuration	Coating
Stiff					
G53211	AS-035145	.035	145	5 cm flexible	PTFE
G53556	AS-038145	.038	145	5 cm flexible	PTFE
G53557	AS-038145-DF	.038	145	5 cm double flexible	PTFE
Extra Stiff					
G53558	AES-025145	.025	145	5 cm flexible	PTFE
G53559	AES-035145	.035	145	5 cm flexible	PTFE
G53560	AES-038145	.038	145	5 cm flexible	PTFE
G53561	AES-035260	.035	260	5 cm flexible	PTFE
G53562	AES-038260	.038	260	5 cm flexible	PTFE
G53563	AES-035145-DF	.035	145	5 cm double flexible	PTFE
G53564	AES-038145-DF	.038	145	5 cm double flexible	PTFE
Ultra Stiff					
G53566	AUS-035145	.035	145	5 cm flexible	PTFE
G53567	AUS-038145	.038	145	5 cm flexible	PTFE
G53568	AUS-035145-DF	.035	145	5 cm double flexible	PTFE
G53569	AUS-038145-DF	.038	145	5 cm double flexible	PTFE



www.cookmedical.com

COOK MEDICAL INCORPORATED

P.O. Box 4195, Bloomington, IN 47402-4195 USA
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll Free Fax: 800.554.8335

COOK IRELAND LTD.

O'Halloran Road, National Technology Park,
Limerick, IRELAND, Phone: +353 61 334440, Fax: +353 61 334441

COOK (CANADA) INC.

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA
Phone: 905.640.7110, Toll Free: 800.668.0300, Fax: 905.640.6804

WILLIAM A. COOK AUSTRALIA PTY. LTD.

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

© COOK 2010 URO-BM-AFCDS-EN-201008

AORTIC
INTERVENTION

CRITICAL
CARE

ENDOSCOPY

INTERVENTIONAL
RADIOLOGY

LEAD
MANAGEMENT

PERIPHERAL
INTERVENTION

SURGERY

UROLOGY

WOMEN'S
HEALTH