

**COOK<sup>®</sup>**  
MEDICAL

# Inserter



## Insert with ease.

### **SCOOP**

Proximal end has a large open trough to allow easy insertion of wire guides and other instrumentation.

### **GRIP**

Exterior grooves provide a convenient grip for additional control and precise insertion.

### **TIP**

Distal tip has a slender taper to ensure a snug fit. Can also be used to puncture the seal on an endoscopic valve. Inserter is compatible with various adapters.



# Inserter



The Cook Medical Inserter has a functional design to facilitate the passage of wire guides and other instrumentation into the working channel of a scope. The Inserter is translucent for visibility and accepts devices up to 3.1 Fr.

Intended for one-time use. Individually packaged sterile. Sold in boxes of 25.

Global Product Number	Order Number	Accepts Instrumentation Fr	Length in
G55179	CI-25	3.1	1.91



[www.cookmedical.com](http://www.cookmedical.com)

**COOK MEDICAL INCORPORATED**

P.O. Box 4195, Bloomington, IN 47402-4195 USA  
Phone: 812.339.2235, Toll Free: 800.457.4500, Toll-Free Fax: 800.554.8335

**COOK (CANADA) INC.**

111 Sandiford Drive, Stouffville, Ontario, L4A 7X5 CANADA  
Phone: 905.640.7110, Toll Free: 800.668.0300

**WILLIAM A. COOK AUSTRALIA PTY. LTD.**

95 Brandl Street, Brisbane Technology Park, Eight Mile Plains  
Brisbane, QLD 4113 AUSTRALIA, Phone: +61 7 3841 1188

**WILLIAM COOK EUROPE ApS**

Sandet 6, DK-4632, Bjaeverskov, DENMARK, Phone: +45 56 86 86 86

© COOK 2010 URO-BM-TIPI-EN-201008

AORTIC INTERVENTION

CRITICAL CARE

ENDOSCOPY

INTERVENTIONAL RADIOLOGY

LEAD MANAGEMENT

PERIPHERAL INTERVENTION

SURGERY

UROLOGY

WOMEN'S HEALTH