

Biopsy efficiently
with a reusable device.



Cook Bx™
BIOPSY DEVICE



www.cookmedical.com

Cook Bx™

BIOPSY DEVICE

The Cook Bx™ Biopsy Device is designed to harvest multiple high-quality histological tissue specimens from the prostate or kidney. The Cook Bx is reusable, which makes it environmentally friendly and minimizes the quantity and cost of disposable waste. The device is designed for use with the Cook Biopsy Needle.

Features

- Easy specimen retrieval that doesn't require removing the needle from the instrument.
- Powerful dual-spring design for rapid-fire action.
- Easy actuation of trigger from front or rear of instrument for single-person operation.
- Automatic safety to prevent accidental firing.
- Small, lightweight design.





Cook Biopsy Needle

The Cook Biopsy Needle is designed for precise placement and sampling. The patented hub design allows for use with the Cook Bx.

- 19 mm sample-notch ensures sufficient tissue sample size.
- Centimeter markings facilitate precise depth placement.
- Echogenic tip enables accurate placement under ultrasound.



Cook Periprostatic Needle

The Cook Periprostatic Needle is used to inject a local anesthetic prior to performing the biopsy.

- Prospective randomized study indicates periprostatic instillation of nerve block is superior to intrarectal instillation.¹
- Can be purchased with or without the biopsy needle for streamlined ordering.

¹ Alavi AS, Soloway MS, Vaidya A, et al. Local anesthesia for ultrasound guided prostate biopsy: a prospective randomized trial comparing 2 methods. *J Urol.* 2001;166(4):1343-1345.

Order Number	Reference Part Number	Needle gage	Length cm	Quantity
Cook Bx Biopsy Device				
G46045	CBX-1			1
Cook Biopsy Needle				
G46044	CBXN-180020	18	20	20
G46071	CBXN-180025	18	25	20
Cook Biopsy Needle and Periprostatic Injection Needle				
G47858	CBXN-180020-PPN	18	20	20
G47859	CBXN-180025-PPN	18	25	20
Periprostatic Injection Needle Only				
G47860	PPN-2020	20	20	20

Not all part numbers shown on this product information sheet may be approved for sale in all regulatory jurisdictions. Consult with your local Cook representative or customer service center for details.

Customer Service Centers - Urology

EMEA: EDI - www.cookmedical.com/edi.do
 Distributors: +353 61239240, ssc.distributors@cookmedical.com
 Austria: +43 179567121, oe.orders@cookmedical.com
 Belgium: +32 27001633, be.orders@cookmedical.com
 Denmark: +45 38487607, da.orders@cookmedical.com
 Finland: +358 972519996, fi.orders@cookmedical.com
 France: +33 171230269, fr.orders@cookmedical.com
 Germany: +49 6950072804, de.orders@cookmedical.com
 Hungary: +36 17779199, hu.orders@cookmedical.com
 Ireland: +353 61239252, ie.orders@cookmedical.com
 Italy: +39 0269682853, it.orders@cookmedical.com
 Netherlands: +31 202013367, nl.orders@cookmedical.com
 Norway: +47 23162968, no.orders@cookmedical.com
 Poland: +48 223060159, pl.orders@cookmedical.com
 Spain: +34 912702691, es.orders@cookmedical.com
 Sweden: +46 858769468, se.orders@cookmedical.com
 Switzerland - French: +41 448009609, fr.orders@cookmedical.com
 Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
 Switzerland - German: +41 448009609, de.orders@cookmedical.com
 United Kingdom: +44 2073654183, uk.orders@cookmedical.com
www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
 E-mail: orders@cookmedical.com

Australia:
 Phone: +61 738411188, 1800777222, Fax: +61 738411288, 1800077283
 E-mail: cau.custserv@cookmedical.com

