

Select from a full line of resectoscope electrodes.



Cook Single-Use RF Electrodes

When you choose a Cook Medical electrode, you're choosing a device that is specifically engineered for the only environment that matters: your patient. The single-use electrodes are indicated for the ablation and coagulation of soft tissue and are intended for use with compatible resectoscopes.



www.cookmedical.com



Cook® Coagulating Perforated Resector

Sold in boxes of 5.



Order Number	Reference Part Number	Fits Fr
Compatible with ACMI Brand USA/Elite Resectoscope		
G46204	PR-2400-AC	24.0
G46247	PR-2600-AC	26.0
G46215	PR-2800-AC	28.0
Compatible with STORZ Brand Single-Stem Resectoscope		
G46221	PR-2400-KS1	24.0
G46222	PR-2700-KS1	27.0
Compatible with STORZ Brand Double-Stem Resectoscope		
G46241	PR-2400-KS2	24.0
G46242	PR-2700-KS2	27.0
Compatible with Olympus Brand OES® Resectoscope		
G46234	PR-2400-OL	24.0
G46235	PR-2600-OL	26.0
G46236	PR-2800-OL	28.0

OES is a registered trademark of Olympus Corporation.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Cook® Collings Knife

Sold in boxes of 5.



Order Number	Reference Part Number	Fits Fr
Compatible with ACMI Brand USA/Elite Resectoscope		
G46208	CK-2400-AC	24.0, 26.0, 28.0
Compatible with STORZ Brand Single-Stem Resectoscope		
G46201	CK-2400-KS1	24.0, 27.0
Compatible with Olympus Brand OES® Resectoscope		
G46227	CK-2400-OL	24.0, 26.0, 28.0

OES is a registered trademark of Olympus Corporation.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Cook® Cutting Loop

Sold in boxes of 5.



Order Number	Reference Part Number	Fits Fr	Wire Diameter inch	Loop Angle degrees
Compatible with ACMI Brand USA/Elite Resectoscope				
G46209	CL-2412-AC	24.0	.012	straight
G46210	CL-2415-AC	24.0	.015	straight
G46211	CL-2612-AC	26.0	.012	straight
G46212	CL-2615-AC	26.0	.015	straight
G46213	CL-2812-AC	28.0	.012	straight
G46214	CL-2815-AC	28.0	.015	straight
Compatible with STORZ Brand Single-Stem Resectoscope				
G46218	CL-2412-KS1	24.0	.012	straight
G46219	CL-2414-KS1	24.0	.014	straight
G46220	CL-2712-KS1	27.0	.012	straight
G46199	CL-2714-KS1	27.0	.014	straight
Compatible with STORZ Brand Double-Stem Resectoscope				
G46203	CL-2414-KS2	24.0	.014	straight
G46240	CL-2714-KS2	27.0	.014	straight
Compatible with Olympus Brand OES® Resectoscope				
G46228	CL-2412-12-OL	24.0	.012	12
G46229	CL-2412-30-OL	24.0	.012	30
G46230	CL-2612-12-OL	26.0	.012	12
G46231	CL-2612-30-OL	26.0	.012	30
G46232	CL-2812-12-OL	28.0	.012	12
G46233	CL-2812-30-OL	28.0	.012	30

OES is a registered trademark of Olympus Corporation.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Cook® Dimpled Vaporization Roller

Sold in boxes of 5.



Order Number	Reference Part Number	Fits Fr	Dimpled Tip Width mm
Compatible with ACMI Brand USA/Elite Resectoscope			
G46205	VR-243M-AC	24.0, 26.0, 28.0	3
Compatible with STORZ Brand Single-Stem Resectoscope			
G46226	VR-243M-KS1	24.0, 27.0	3
Compatible with STORZ Brand Double-Stem Resectoscope			
G46246	VR-243M-KS2	24.0, 27.0	3
Compatible with Olympus Brand OES® Resectoscope			
G46239	VR-243M-OL	24.0, 26.0, 28.0	3


OES is a registered trademark of Olympus Corporation.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Cook® Roller Ball and Roller Barrel


Sold in boxes of 5.



Order Number	Reference Part Number	Fits Fr	Tip Configuration
Compatible with ACMI Brand USA/Elite Resectoscope			
G46216	RB-243M-AC	24.0, 26.0, 28.0	3 mm ball
G46217	RB-245M-AC	24.0, 26.0, 28.0	5 mm ball
G46202	RL-243M-AC	24.0, 26.0, 28.0	3 mm barrel
Compatible with STORZ Brand Single-Stem Resectoscope			
G46223	RB-243M-KS1	24.0, 27.0	3 mm ball
G46224	RB-245M-KS1	24.0, 27.0	5 mm ball
G46225	RL-243M-KS1	24.0, 27.0	3 mm barrel
Compatible with STORZ Brand Double-Stem Resectoscope			
G46243	RB-243M-KS2	24.0, 27.0	3 mm ball
G46244	RB-245M-KS2	24.0, 27.0	5 mm ball
G46245	RL-243M-KS2	24.0, 27.0	3 mm barrel
Compatible with Olympus Brand OES® Resectoscope			
G46200	RB-243M-OL	24.0, 26.0, 28.0	3 mm ball
G46237	RB-245M-OL	24.0, 26.0, 28.0	5 mm ball
G46238	RL-243M-OL	24.0, 26.0, 28.0	3 mm barrel

OES is a registered trademark of Olympus Corporation.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.



Customer Service

EMEA: EDI - www.cookmedical.com/edi.do
Distributors: +353 61239240, ssc.distributors@cookmedical.com
Austria: +43 179567121, oe.orders@cookmedical.com
Belgium: +32 27001633, be.orders@cookmedical.com
Denmark: +45 38487607, da.orders@cookmedical.com
Finland: +358 972519996, fi.orders@cookmedical.com
France: +33 171230269, fr.orders@cookmedical.com
Germany: +49 6950072804, de.orders@cookmedical.com
Hungary: +36 17779199, hu.orders@cookmedical.com
Ireland: +353 61239252, ie.orders@cookmedical.com
Italy: +39 0269682853, it.orders@cookmedical.com
Netherlands: +31 202013367, nl.orders@cookmedical.com
Poland: +48 223060159, pl.orders@cookmedical.com
Spain: +34 912702691, es.orders@cookmedical.com
Sweden: +46 858769468, se.orders@cookmedical.com
Switzerland - French: +41 448009609, fr.orders@cookmedical.com
Switzerland - Italian: +41 448009609, it.orders@cookmedical.com
Switzerland - German: +41 448009609, de.orders@cookmedical.com
United Kingdom: +44 2073654183, uk.orders@cookmedical.com

www.cookmedical.com

Americas: EDI - www.cookmedical.com/edi.do
Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335
E-mail: orders@cookmedical.com

Australia:
Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283
E-mail: cau.custserv@cookmedical.com

