



# Prepare for suprapubic drainage.

**Universa<sup>®</sup>**  
MALECOT DRAINAGE CATHETER

The Universa Malecot Drainage Catheter Suprapubic Set is used for percutaneous urine drainage from the genitourinary system. Universa products are crafted with quality materials to give you the fit and feel you're accustomed to and to meet your routine, everyday drainage needs.




# Universa®

## MALECOT DRAINAGE CATHETER

The set includes:

- Malecot catheter
- Locking needle
- Connecting tube with one-way stopcock\*



| Order Number                     | Reference Part Number | Catheter Fr | Catheter Length cm | Needle gage |
|----------------------------------|-----------------------|-------------|--------------------|-------------|
| <b>Malecot Drainage Catheter</b> |                       |             |                    |             |
| Suprapubic Set                   |                       |             |                    |             |
| G56334                           | UMCP-8-15             | 8.0         | 15                 | 19          |
| G55821                           | UMCP-8-30             | 8.0         | 30                 | 19          |
| G56335                           | UMCP-10-15            | 10.0        | 15                 | 18          |
| G55822                           | UMCP-10-30            | 10.0        | 30                 | 18          |
| G55823                           | UMCP-12-30            | 12.0        | 30                 | 16          |
| G55824                           | UMCP-14-30            | 14.0        | 30                 | 15          |
| G55825                           | UMCP-16-30            | 16.0        | 30                 | 15          |

\*The one-way stopcock is not available on 16 Fr catheters.

Some products or part numbers may not be available in all markets. Contact your local Cook representative or Customer Service for details.

### Customer Service

EU Website: [cookmedical.eu](http://cookmedical.eu)  
 EDI: [cookmedical.eu/edi](http://cookmedical.eu/edi)  
 Distributors: +353 61239240, [ssc.distributors@cookmedical.com](mailto:ssc.distributors@cookmedical.com)  
 Austria: +43 179567121, [oe.orders@cookmedical.com](mailto:oe.orders@cookmedical.com)  
 Belgium: +32 27001633, [be.orders@cookmedical.com](mailto:be.orders@cookmedical.com)  
 Denmark: +45 38487607, [da.orders@cookmedical.com](mailto:da.orders@cookmedical.com)  
 Finland: +358 972519996, [fi.orders@cookmedical.com](mailto:fi.orders@cookmedical.com)  
 France: +33 171230269, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Germany: +49 6950072804, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 Ireland: +353 61239252, [ie.orders@cookmedical.com](mailto:ie.orders@cookmedical.com)  
 Italy: +39 0269682853, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Netherlands: +31 202013367, [nl.orders@cookmedical.com](mailto:nl.orders@cookmedical.com)  
 Poland: +48 223060159, [pl.orders@cookmedical.com](mailto:pl.orders@cookmedical.com)  
 Spain: +34 912702691, [es.orders@cookmedical.com](mailto:es.orders@cookmedical.com)  
 Sweden: +46 858769468, [se.orders@cookmedical.com](mailto:se.orders@cookmedical.com)  
 Switzerland - French: +41 448009609, [fr.orders@cookmedical.com](mailto:fr.orders@cookmedical.com)  
 Switzerland - Italian: +41 448009609, [it.orders@cookmedical.com](mailto:it.orders@cookmedical.com)  
 Switzerland - German: +41 448009609, [de.orders@cookmedical.com](mailto:de.orders@cookmedical.com)  
 United Kingdom: +44 2073654183, [uk.orders@cookmedical.com](mailto:uk.orders@cookmedical.com)

USA Website: [cookmedical.com](http://cookmedical.com)  
 EDI: [cookmedical.com/edi.do](http://cookmedical.com/edi.do)  
**Americas:**  
 Phone: +1 812.339.2235, 800.457.4500, Fax: 800.554.8335  
 E-mail: [orders@cookmedical.com](mailto:orders@cookmedical.com)  
**Australia:**  
 Phone: +61 734346000, 1800777222, Fax: +61 734346001, 1800077283  
 E-mail: [cau.custserv@cookmedical.com](mailto:cau.custserv@cookmedical.com)



CC, URO-85X11