

G57703 is shown here.
Product components may vary
by part number.

Blue Rhino® G2-Multi

PERCUTANEOUS TRACHEOSTOMY INTRODUCER TRAY

Tray components

- ① Blue Rhino G2-Multi Percutaneous Tracheostomy Dilator
 - ② Tracheostomy tube loading dilators¹
 - ③ 15 gage introducer needle
 - ④ 15 gage introducer needle with FEP sheath
 - ⑤ Safe-T-J® Fixed Core Wire Guide
 - ⑥ Introducer dilator
 - ⑦ Guiding catheter
 - ⑧ 10 4 x 4 inch gauze sponges
 - ⑨ 6 and 12 mL Monoject® syringes
 - ⑩ Needle holder cup
 - ⑪ Sterile lubricating jelly²
 - ⑫ Wire guide dispenser with J-tip straightener
 - ⑬ #15 safety scalpel
 - ⑭ 2 5mL ampoules of 1.5% lidocaine HCL²
 - ⑮ 22 and 25 gage Monoject needles
 - ⑯ Double swivel connector
 - ⑰ 30 x 30 inch CSR wrap
 - ⑱ 4 x 4 inch drain sponge
 - ⑲ 5 inch long curved hemostat
 - ⑳ 5 inch long serrated hemostat
 - ㉑ Tracheostomy tube holder
 - ㉒ 44 x 30 inch folded drape
 - ㉓ 2/0 polypropylene blue monofilament suture
 - ㉔ Two Chloraprep® with tint applicators²
 - ㉕ Filter straw²
- Tracheostomy tube with disposable inner cannula (not pictured)³

1. Size and quantity of tracheostomy tube loading dilators will vary based on tray or set configurations.
2. Not included in trays without pharmaceuticals.
3. Included only where listed.



Monoject is a registered trademark of Medtronic.

Chloraprep is a registered trademark of CareFusion 2200, Inc.

